

Notes: _____

Types of Provox HME Attachments

Provox Adhesives

- Are available in 2 shapes and 4 adhesive types.
- Are taped onto the skin around your stoma.



Provox Adhesives

Provox LaryTube

- Helps to prevent the stoma from shrinking.
- Can be secured to the neck with a Provox TubeHolder, Provox LaryClips, or a Provox Adhesive.
- Is available in 3 different types with various lengths and widths.



Provox LaryTube Fenestrated with Provox TubeHolder

Provox LaryButton

- Is a soft stoma button that fits comfortably in the stoma.
- Can be secured to the neck with the Provox TubeHolder or Provox LaryClips.
- Is available in various lengths and widths.



Provox LaryTube with Ring in a Provox XtraBase adhesive

Tip: Prior to choosing an adhesive or other attachment, it is important to meet with your clinician who can assist you with examining your stoma and selecting the appropriate attachment method.



Provox LaryButton with Provox TubeHolder



Provox LaryButton with Provox LaryClips and Provox FreeHands HME